



www.free-foundation.org

shower chair

A qualified medical professional should educate and train users for appropriate use. It is your responsibility to use this device correctly.

A shower chair can be used not only by the injured or disabled person in taking a shower, but also as a comfortable device for everyone. So, you must check the chair regularly, and its stability to prevent accidents.

Installation

A qualified medical professional should install the shower chair to ensure appropriate installation. Place the chair in the tub or shower stall. Be sure it's stable. Remove obstructions that can topple and hurt someone or cause them to slip. Sit in the center of the chair, never on the edge. Most of the shower chairs have adjustable legs. You can adjust the height so users can sit and stand up easily. Make sure that the chair stands on a level surface and with all legs adjusted to the same height. The rubber tips on your shower chair should be inspected regularly. Worn or damaged tips should be replaced immediately. The security of all adjustment mechanisms should be checked frequently.

It is your responsibility to maintain and provide upkeep of this item. It is your responsibility to use this device correctly.