



www.free-foundation.org



clamp on grab bar

A qualified medical professional should educate and train users for appropriate use. It is your responsibility to use this device correctly.

Providing "something to hold onto" and provides stability. Grab bars cannot be safely installed on fiberglass or plastic bathtubs. It is important to select a grab bar with the correct type of mounting bracket to provide a tight and proper fit.

The security of the mounting should be tested by pulling forcefully on the bar in all directions. The mounting of any grab bar should be re-checked frequently for tightness and security. A loose or unstable grab bar is worse than no grab bar at all. Contact a qualified medical professional to secure any loose or unstable bar. Keeping the grab bar clean and free of soap film, oil, or grease will help to retain any no-slip characteristic that may have been provided by the manufacturer.

It is your responsibility to maintain and provide upkeep of this item. It is your responsibility to use this device correctly.