



[www.free-foundation.org](http://www.free-foundation.org)



**geri chair**

**A qualified medical professional should educate and train users for appropriate use. It is your responsibility to use this device correctly.**

Geri chairs are specialized recliners that are upholstered in non-permeable, easily sanitized vinyl. This makes the chair easy to clean and provides firm support for users. Additionally, geri chairs have wheels on their feet, which makes transporting individuals easy. Many chairs come with positioning cushions and activity tray tables.

These chairs are a perfect addition in the home or in a clinical environment where some individuals need more support than a traditional wheelchair can provide. To recline, most Geri chairs have a bar or lever you use to move the chair in the desired position and lock it in place. Similarly, the chairs should have casters with locks to easily move and secure the recliner. These casters should always be locked when transferring in or out of the chair or when stationary.

**It is your responsibility to maintain and provide upkeep of this item. It is your responsibility to use this device correctly.**